


Durable by design: Zenith Alpha[®] Abdominal

Learn more about the details that went into the device design.

Proven design.

From the proximal to the distal seal, the Zenith Alpha Abdominal device is engineered to help you deliver a durable repair.



Intuitive introduction system.

Zenith Alpha Abdominal features a new ergonomic handle that simplifies device deployment while maintaining precision and control.



1 Unsheathe



2 Unlock




3 Rotate




Patient-specific fit.

With a broad range of graft lengths and diameters, and a low profile introduction system, Zenith Alpha Abdominal is designed to help you fit the device to the patient.



The low profile introduction system helps improve vessel access.



The hydrophilic-coated Flexor[®] Introducer Sheath helps the device track through iliac anatomy.

