by part number. 20 POLYPROPYLENE 23 TOSTOUCH 3 COVIDIEN" cilon" monoject REF 10.2 cm x 10.2 cm) es absorbantes COOK

Blue Rhino G2-Multi

PERCUTANEOUS TRACHEOSTOMY INTRODUCER TRAY

Tray components

- 1 Blue Rhino G2-Multi Percutaneous Tracheostomy Dilator
- 2 Tracheostomy tube loading dilators¹
- 3 15 gage introducer needle
- 4 15 gage introducer needle with FEP sheath
- 5) Safe-T-J® Fixed Core Wire Guide
- (6) Introducer dilator
- 7 Guiding catheter
- 8 10 4 x 4 inch gauze sponges
- 9 6 and 12 mL Monoject® syringes
- 10 Needle holder cup
- Wire guide dispenser with J-tip straightener
- (12) #15 safety scalpel
- 13) 22 and 25 gage Monoject needles
- 14) Double swivel connector
- (15) 30 x 30 inch CSR wrap
- (16) 4 x 4 inch drain sponge
- 5 inch long curved haemostat
- (18) 5 inch long serrated haemostat
- 19 Tracheostomy tube holder
- 20 44 x 30 inch folded drape
- 2/0 polypropylene blue monofilament suture

Tracheostomy tube with disposable inner cannula (not pictured)²

- 1. Size and quantity of tracheostomy tube loading dilators will vary based on tray or set configurations.
- 2. Included only where listed.

