

Blue Rhino[®] G2-Multi

PERCUTANEOUS TRACHEOSTOMY INTRODUCER SET

Set components

- 1 Blue Rhino G2-Multi Percutaneous Tracheostomy Dilator
- 2 Tracheostomy tube loading dilators¹
- (3) 15 gage introducer needle
- 4 15 gage introducer needle with FEP sheath
- (5) Safe-T-J[®] Fixed Core Wire Guide
- 6 Introducer dilator
- 7 Guiding catheter
- 8 Four 4 x 4 inch gauze sponges
- 9 6 mL Monoject® syringe
- (10) Needle holder cup
- (11) Sterile lubricating jelly
- (12) Wire guide dispenser with J-tip straightener
- (13) #15 safety scalpel
- (14) 5 inch long curved haemostat

Tracheostomy tube with disposable inner cannula (not pictured)²

 Size and quantity of tracheostomy tube loading dilators will vary based on tray or set configurations.
Included only where listed.

